

PROEIZ™

MESO RF

Excellent anti-aging treatment on the eye area and face!

- Facial Lifting
- Skin Tightening
- Wrinkles Reduction
- Pores Reduction
- Skin Rejuvenation
- Dark Circles
- Eye Contour Lifting
- Eye Bags & Edema

Eye problems occur at all times

Work pressure, lack of sleep, environmental pollution, and others caused skin problems

Dark Circles



Before

After


Eye Bags & Edema



Before

After

Wrinkles



Before

After






Meso RF Technology




1. Nutrients Permeation

Deep permeation function with professional skin care products, this is to stimulate the skin, promote collagen regeneration and soft tissue reconstruction, beautify skin complexion. An excellent anti-aging effect!




2. High Frequency Vibration

Further accelerate the blood circulation, promote metabolism.



3. Red Light Care

630 nm red light, with a strong penetration, can penetrate into skin tissue, repair cells, and activate biological macromolecule activity.




4. Radio Frequency

Deep skin heating, stimulate collagen regeneration, restore skin elasticity and to achieve long-term purpose of reducing skin wrinkles.


Skin care products effect amplifier, high nutrients absorption rate up to 98%

Traditional Way of Permeation



Skin care products mainly in the superficial skin layer, nutrient absorption rate of less than 15%

NEW Smart Transdermal System



With Meso RF, skin care products break through the skin, go direct to dermis, improve nutrient absorption up to 98%

Meso RF Treatment Results



Before Treatment



After Treatment
(4 weeks, once a week)


Before Treatment:

There are eye bags under the eyes, tired looking with dark trench.


After Treatment:

After the first use, the skin is **highly moisturised**. In a week, the bags under the eyes are significantly smaller. After a month, the face is firmer, and bags under the eyes are obviously faded and **the effect is really good!**


Smart Transdermal System



1 Rotate the head to open



2 Remove the solution bottle



3 Open tip of bottle



4 Pour in treatment solution

Technical Specifications

| | | |
|-----------------------|---|---|
| Mode | : | RF, Vibration, Nutrient Permeation, Red Light 630nm |
| Vibration Level Range | : | Level 1 - 100Hz, Level 2 - 125Hz, Level 3 - 150Hz |
| RF Level Range | : | Level 1 - 15W, Level 2 - 20W, Level 3 - 25W |
| Strength | : | 0.1W |
| Power Supply | : | AC (100-240)V, (50-60)Hz |
| Input Power | : | Maximum 30W |
| Output Power | : | Maximum 25W |
| Net Weight | : | 1.4kg |
| Total Weight | : | 2.4kg |
| Dimensions | : | Length 25.3cm x Width 25.3cm x Height 16cm |
| Packing Dimensions | : | Length 35.5cm x Width 31.5cm x Height 16cm |

PROEIZ™

MESO RF

优越的抗衰老治疗眼部
和脸部！

- ☑ 面部提升
- ☑ 紧致皮肤
- ☑ 减少皱纹
- ☑ 减少毛孔
- ☑ 嫩肤
- ☑ 黑眼圈
- ☑ 眼部提升
- ☑ 眼袋和水肿

每时每刻眼睛都能出现问题

工作压力，睡眠不足，环境污染等
导致皮肤问题

黑眼圈



之前



之后

眼袋和水肿



之前



之后

皱纹



之前



之后






Meso RF 美容仪技术原理




1. 营养液在线导入

深层导入功能，配合专业护肤品，激发皮肤，深层胶原再生，软组织重建肌肤得以新生，达到美肤，抗衰效果！



2. 高频振动

进一步加速血液循环，促进新陈代谢。



3. 红光波护理

630 纳米红光，具有很强的穿透性，可以渗透到人体组织中，修复人体细胞，激活生物大分子活性。




4. 射频

皮肤真皮加热，刺激胶原再生，使皮肤弹性组件恢复，达到长期减少皮肤皱纹的目的。

护肤品效果放大器，营养吸收程度高达 **98%**


传统导入方式



护肤品主要在浅表层滋养，营养吸收不及 **15%**



智能在线导入



Meso RF 美容仪，使护肤品突破表皮，直达真皮层，营养吸收程度高达 **98%**

Meso RF 美容仪治疗效果



治疗前

治疗后
(4周, 1周1次)


治疗前:

有眼袋, 疲惫的暗沉泪沟。


治疗后:

第一次用后皮肤就特别滋润。一个星期眼袋就明显变小了。一个月后脸上的皮肤紧实了, 眼袋泪沟也明显变淡了, 这效果真的不错!


护肤品智能在线导入



1 沿箭头旋转头部



2 将头部从手柄内拔出



3 将手柄头拔出



4 注入护肤品

技术规格

| | | |
|-----------------------|---|---|
| Mode | : | RF, Vibration, Nutrient Permeation, Red Light 630nm |
| Vibration Level Range | : | Level 1 - 100Hz, Level 2 - 125Hz, Level 3 - 150Hz |
| RF Level Range | : | Level 1 - 15W, Level 2 - 20W, Level 3 - 25W |
| Strength | : | 0.1W |
| Power Supply | : | AC (100-240)V, (50-60)Hz |
| Input Power | : | Maximum 30W |
| Output Power | : | Maximum 25W |
| Net Weight | : | 1.4kg |
| Total Weight | : | 2.4kg |
| Dimensions | : | Length 25.3cm x Width 25.3cm x Height 16cm |
| Packing Dimensions | : | Length 35.5cm x Width 31.5cm x Height 16cm |