

## Body Slimming & Shaping

- US Cavitation - Fat Tissue Reduction**  
 Cavitation Principle causes thousands of microscopic air bubbles to surround the fat cell membranes with resonant frequency, causing enough pressure for fat cells to break up instantly.
- OCTA-POLE RF + Laser Lipolysis & TEN-POLE RF - Sculpting Figure**  
 Radio Frequency can tighten muscles of the arms, legs, thighs, buttocks, back, abdominal muscles, re-shaping body contour.
- Vacuum + RF - Stimulate Metabolism, Lymphatic Drainage**  
 Vacuum Principle can improve the blood circulation, do lymphatic drainage, promote tissue metabolism.
- Cryo Therapy - Shape, Firm & Tighten Skin**  
 Cryo Therapy reduces fat cells and tighten skin.

## Face Lift & Eye Contour Lift

- TRIPOLAR RF (Face) & BIPOLAR RF (Eye) - Restores Elasticity, Firms Skin**  
 Radio Frequency energy penetrates the skin and stimulates the contraction of collagen and the production of new collagen for skin tightening.



US Cavitation

OCTA-POLE RF + Laser

TEN-POLE RF

Vacuum + RF

Cryo Therapy

TRIPOLAR RF

BIPOLAR RF



### Technical Specifications

Voltage	: 220V - 240V
Output Power	: 500W
Screen	: 7"inch TFT Touch Screen
Cavitation Frequency	: 40KHz
Cryo Temperature	: 5°C
Laser Wave Length	: 650nm
RF Frequency	: 5MHz
Working Temperature	: 10°C - 40°C

## 纤体塑身

- US Cavitation - 脂肪组织减少**  
 气蚀原理导致数以千计的微小气泡与谐振频率包围脂肪细胞膜，造成足够的压力使脂肪细胞瞬间分裂。
- OCTA-POLE RF + Laser Lipolysis & TEN-POLE RF - 雕塑身型**  
 无线电频率可以收紧手臂，小腿，大腿，臀部肌肉，背部，腹部肌肉，重新塑造身体轮廓。
- Vacuum + RF - 促进循环和新陈代谢**  
 真空原理可改善血液循环，做淋巴排毒，促进组织的新陈代谢。
- Cryo Therapy - 塑形，结实及紧致肌肤**  
 冷冻疗法减少脂肪细胞和紧实肌肤。

## 面部与眼部提升

- TRIPOLAR RF (面部) & BIPOLAR RF (眼部) - 恢复弹性，紧致肌肤**  
 射频能量穿透皮肤，刺激胶原蛋白的收缩和生产新胶原蛋白，用于皮肤收紧。



US Cavitation

OCTA-POLE RF + Laser

TEN-POLE RF

Vacuum + RF

Cryo Therapy

TRIPOLAR RF

BIPOLAR RF

### 技术指标

电压	: 220V - 240V
电源输出	: 500W
屏幕	: 7" TFT 触摸屏
气蚀频率	: 40KHz
冷冻疗法温度	: 5°C
激光波长	: 650nm
无线电频率	: 5MHz
使用温度	: 10°C - 40°C

